

<u>Set-Up</u>

- Shape the stylet with a slight bend at the proximal cuff. Bend angle is typically 35'.
- To prevent fogging, apply anti-fog or warm the stylet by using a warm blanket or immersing in a warm saline bottle.
- Cut ET tube so that the stylet is recessed inside the tube. Place ET tube into the tube stop.
- Set the oxygen flow through the O2 port.

ACTIONS	SIDE VIEW	PLACEMENT	USER'S VIEW
 Perform direct laryngoscopy exposing epiglottis edge Option: If glottic opening is visible Intubate under direct vision using Levitan as normal stylet 	THE REAL		
 Using direct vision, place the stylet in position Place the tip of the stylet under the epiglottis and away from mucosa. Do not advance stylet's tip beyond your direct vision 		and and a	
 Switch from Direct view to fiberoptic view without moving stylet View the epiglottis edge, posterior cartilage and larynx Gentle maneuvers (rocking back) may assist in the event that the cords are not immediately in view 	The second second	1 A	
 Slowly advance the scope under fiberoptic visualization through the vocal cords Remove the laryngoscope Slide the ET tube off with left hand Visualize the tube entering the trachea 		A A	

NOTE: The Levitan is versatile and can be used many ways. See how-to guide for alternative techniques on back.



13355 10th Ave N. Minneapolis, MN 55441 (800) 359 - 2372 Doc# 910091-001-A

Alternative Techniques

WITHOUT A LARYNGOSCOPE

Useful for patients with: Routine & Difficult Intubations, Cervical Spine & Dental Concerns, Limited Mouth Opening

Change the curvature of stylet as desired
 Create an open channel

 Lift tongue with 4x4 of cotton
 Jaw lift/Jaw thrust

2 Insert the stylet midline (or retro-molar approach) do not advance beyond base of tongue
Start visualizing early
Pass the tongue, epiglottis, and finally pass the stylet through the cords

•While maintaining visualization, withdraw the stylet and advance the ET tube

e ET

CRICOTHYROTOMY

- This technique provides
- •Visual confirmation
- •Optically guided training of cricothyrotomy procedure



1: Paladino L, DuCanto J, Manoach S. Development of a rapid, safe, fiber-optic guided, single-incision cricothyrotomy using a large ovine model: a pilot study. Resuscitation. 2009 Sep;80(9):1066-9. Epub 2009 Jul 15.

ASSISTED BY ORAL AIRWAY

Oral airways assist Levitan intubations by creating a more open and navigable channel

"If you can't see the vocal cords after getting the scope through the airway, move the airway with your left hand as you look through the eyepiece and the laryngeal inlet will appear:"

- Eric Flaten, Ridgecrest, CA



WITH A VIDEO LARYNGOSCOPE

The Levitan is used as an adjunct to steer the tube around corners and provides visual confirmation of tube placement

•Anterior Airways

•Cervical Spine & Dental concerns - Use less force with the laryngoscope while the Levitan Stylet's curvature provides the view to the cords



THROUGH LARYNGEAL MASK AIRWAY

This technique allows for:

•Visual confirmation of LMA/SGA and ETT placement •Use the stylet's rigidity to maneuver LMA/SGA •Isolate mucus / vomit during an fiberoptic intubation •Option to wake patient with LMA/SGA

How To:

 1. Lubricate ETT, Scope and interior of the ILMA/SGA

 Albert Einstein MediatemetrLevitan Stylet Distal Tip 40 degrees

- 3. Remove (and keep) the connector from the ILMA and insert the ETT-Levitan combination through the ILMA
 - 4. Look into Levitan and visualize
 - Pass through & visualize mask of ILMA SGA
 Visualize vocal cords and pilot the EVE / Levitan through the vocal cords
 - 5. Gently advance the ETT off the stylet to complete the intubation





A Wonderfully Simple Guide Sheet To The Use Of The Clarus Levitan Stylet With The Cook Gas ILA, Version 2.0. Intubation Through the Cook ILA With and Without an Optical Stylet

Intubation Through the Cook ILA With and Without an Optical Stylet Richard M. Levitan1, William C. Kinkle2, Kenneth Butler3, James S.