

# CLARUS<sup>®</sup>

## PEDIATRIC SHIKANI



Malleable stylet can be shaped  
to fit your patient's difficult airway

When time and safety are critical, the Clarus Pediatric stylet is a portable and durable tool for a swift and conclusive intubation.

- Sized specifically for delicate pediatric airways
- Fiberoptic technology provides superior vision for a gentle intubation
- Reusable and versatile, giving you incomparable value
- Works in conjunction with many LMA's
- Portable and reusable stylet that shapes to fit your patient's difficult airway

SUPERIOR DESIGN, SUPERIOR ACCESS

 **CLARUS<sup>®</sup>**  
MEDICAL

# CLARUS PEDIATRIC SHIKANI

30001-10



The Clarus Pediatric is designed to give you superior access to the pediatric airway. With its malleable stylet and atraumatic tip, the Clarus Pediatric gives you the size and features needed for a successful intubation. Optimum grade fiberoptics and a high-resolution eyepiece make this stylet a must for those working with delicate airways.

- Conveniently sized for portability
- ① High-resolution eyepiece for confident placements
- ② Port for oxygen insufflation
- ③ Malleable stylet accommodates ET tubes size 5.5 and below
- ④ Atraumatic tip helps protect sensitive airways
- ⑤ Soakable stylet and light source allow for safe and efficient sterilization

1. Shukry, M. et al. Case Report: Management of the difficult airway with Shikani Optical Stylet. *Pediatric Anesthesia*. 2005; 15: 342-345.
2. Stricker, P., Fiadjoe, J.E., McGinnis, S. Intubation of an infant with Pierre Robin sequence under dexmedetomidine sedation using the Shikani Optical Stylet. *Acta Anesthesiologica Scand*. 2008 Jul; 52(6): 866-867.
3. Pfitzner, L., Cooper, M.G., Ho, D. The Shikani Seeing Stylet for difficult intubation in children: Initial experience. *Anaesth. Intensive Care*. 2002; 30: 462-466.

Stylet Length	26.9cm	<b>Cleaning:</b> The Shikani stylet and Turbo LED are fully submersible. Common cleaning procedures include Steris, Sterrad or Cidex (No Autoclave). See full Directions for Use for additional cleaning instructions.
Stylet Diameter	2.4mm	
ET Tube Size Range	2.5 - 5.0	
O2 Port Insufflation	As Needed	

[www.clarus-medical.com/airway](http://www.clarus-medical.com/airway)

**CLARUS**<sup>®</sup>  
MEDICAL

13355 10th Ave N. Suite 110  
Minneapolis, MN 55441 USA  
airway@clarus-medical.com

763-525-8450

800-359-2372

910073-004-A  
US Patents: 5735792, 5512034, 5941816